

DiaExpert Sensor

continuous glucose
monitoring system

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ENGLISH

Quick Start Guide

Get started with Continuous Glucose Monitoring (CGM) quickly and easily with this Quick Start Guide.

READ THIS INSERT AND ALL OF THE LABELLING PROVIDED WITH THE CGM APP BEFORE HANDLING THE SENSOR KIT.

Product name: Continuous Glucose Monitoring System Sensor
Product Model: GX-01E
REF: C046000, C046001

Indication for use

The CGM sensor is a real time, continuous glucose monitoring device. When the system is used together with compatible devices, it is indicated for the management of diabetes in adults (age 18 and older). It is designed to replace fingerstick blood glucose testing for diabetes treatment decisions. Interpretation of the system results should be based on the glucose trends and several sequential readings over time. The system also detects trends and tracks patterns, and aids in the detection of episodes of hyperglycemia and hypoglycemia, facilitating both acute and long-term therapy adjustment.

STEP 1: INSTALL DIAEXPERT APP

Search DiaExpert on the Google Play or App store.

STEP 2: REGISTER AN ACCOUNT

Use an e-mail account to register a DiaExpert account.

STEP 3: APPLY THE SENSOR

The sensor is located inside the sensor applicator. Follow the instructions to prepare and apply the sensor on the back of your upper arm. The sensor has a small, flexible tip that is inserted just under the skin. The sensor can be worn for up to 15 days.

Description

1. Select the back of the upper arm as the insertion site. (avoid the muscle on the outer side of the upper arm.)
2. Sterilize the insertion site with an alcohol wipe and let it dry completely.
3. Hold the applicator vertically with the cover facing up. Unscrew the cover of the applicator and make sure to take out the desiccant block before inserting the sensor.
4. Align the opening of the applicator with the insertion site tightly and press the white implantation button of the applicator. After hearing the sound of the spring retracting, the sensor will be inserted underneath the skin and the puncture needle will automatically retreat back into the applicator.
5. Gently pull the sensor applicator away from the body, and the sensor should now be attached to the skin.
6. After installing the sensor, make sure that the sensor is firmly in place. Use your fingers to flatten the edge of the sensor to avoid wrinkles and warping.

STEP 4: PAIR SENSOR WITH DIAEXPERT APP

1. Log in on the DiaExpert app.
 2. Make sure the phone's bluetooth is turned on.
 3. Click "pairing" on the home dashboard.
 4. Select and click the correct SN number of the sensor listed on the Nearby Device list. Make sure the selected SN number is in accordance with the number printed on the packaging box and the applicator cover.
 5. Wait for a successful pairing.
- Refer to the User Guide for more details of pairing procedures and the subsequent use. The electronic manual should be accessible for viewing using a smartphone or computer with a web browser and PDF viewing capability. If you choose to view the electronic manual on the DiaExpert App, make sure to download the application on the compatible device models listed in our equipment list. Please refer to the following website to check the compatible device list and User Guide:

<https://www.diaexpert.com>

Contraindications

The CGM sensor must be removed prior to Magnetic Resonance Imaging (MRI). Don't wear your CGM sensor for computed tomography (CT) scan, or high-frequency electrical heat (diathermy) treatment. Taking higher than the maximum dose of acetaminophen (e.g. > 1 gram every 6 hours in adults) may affect the CGM readings and make them look higher than they really are. The CGM system was not evaluated for the following persons:

- Pregnant women
- Peritoneal dialysis patients
- Patients with implanted pacemakers
- Patients with coagulation disorders or those taking anticoagulant drugs

Precautions

- No modifications to the CGM system are allowed. Unauthorized modification of the CGM may cause the product to malfunction and become unusable.
- Before using this product, you need to read the Instruction Manual or be trained by a professional. No doctor's prescription is required for use at home.
- The CGMS contains many small parts that can be dangerous if swallowed.
- During rapid changes in blood glucose (more than 0.1 mmol/L or 1.8mg/dL per minute), glucose levels measured in interstitial fluid by the CGMS may not be the same as blood glucose levels. When blood glucose levels drop rapidly, the sensor may produce a higher reading than the blood glucose level; Conversely, when blood glucose levels rise rapidly, the sensor may produce a lower reading than the blood glucose level. In these cases, the sensor's reading is checked by a fingertip blood test using a glucose meter.
- When it is necessary to confirm hypoglycemia or near-hypoglycemia as measured by a glucose sensor, a fingertip blood test should be performed using a glucose meter.
- Severe dehydration or excessive loss of water may result in inaccurate results. When you suspect you are dehydrated, consult a healthcare professional immediately.
- If you think the CGMS sensor reading is inaccurate or inconsistent with the symptoms, use a blood glucose meter to test your blood glucose level or calibrate the glucose sensor. If the problem persists, remove and replace the sensor.
- The performance of the CGMS has not been evaluated when used with another implantable medical device, such as a pacemaker.
- Details of what interferences may affect the accuracy of the detection are given in "Potential Interference information".
- The sensor loosens or takes off may cause the APP to have no readings.
- If a sensor tip breaks, do not handle it yourself. Please seek professional medical help.
- This product is waterproof and can be worn during showers and swimming, but do not bring sensors into the water more than 2 meters deep for longer than 1 hour.
- CGMS readings should only be used as a reference for the supplemental monitoring of diabetes mellitus and should not be used as a basis for clinical diagnosis.
- While extensive user testing was done on Dia Expert CGMS in Type 1 and Type 2 diabetic patients, the study groups did not include women with gestational diabetes.
- If the product is not working properly or has been damaged, stop using the product.
- For user safety, storage, disposal and handling, please refer to the system instruction for use.

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