

DiaExpert Sensor 2

Compatible smartwatch list

The CGM 2 App can also pair with a supported Apple Watch to display glucose data on your wrist. Please review the minimum system requirements below, or confirm your model against the compatible list.

› Minimum Requirements for Software Installation

Requirement	Apple Watch
Operating System (OS)	Watch OS 10 and above
Memory	1GB RAM
Storage	Minimum 200 MB
Network	WLAN or cellular network, as well as Bluetooth function
Screen Resolution	324 x 394 pixels and above

› Compatible Apple Watch Models

This list reflects Apple’s official watchOS compatibility and may be updated as new models are released.

Apple Watch Series 4	Apple Watch SE (1st Gen)	Apple Watch Ultra
Apple Watch Series 5	Apple Watch SE (2nd Gen)	Apple Watch Ultra 2
Apple Watch Series 6	Apple Watch SE 3	Apple Watch Ultra 3
Apple Watch Series 7		
Apple Watch Series 8		
Apple Watch Series 9		
Apple Watch Series 10		
Apple Watch Series 11		

Note: A smartwatch is a companion accessory. Even when a supported Apple Watch is used, the CGM 2 App must still be installed on a compatible iPhone to receive and display sensor data.